

DGKQ Rabbit Polyclonal Antibody

Catalog No: #54797

Package Size: #54797-1 50ul #54797-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DGKQ Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human DGKQ (NP_001338.2).
Other Names	DGKQ;DAGK;DAGK4;DAGK7
Accession No.	Uniprot:P52824GeneID:1609
Calculated MW	101kDa
SDS-PAGE MW	101kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DGKQ antibody.

Background

The protein encoded by this gene contains three cysteine-rich domains, a proline-rich region, and a pleckstrin homology domain with an overlapping Ras-associating domain. It is localized in the speckle domains of the nucleus, and mediates the regeneration of phosphatidylinositol (PI) from diacylglycerol in the PI-cycle during cell signal transduction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.