

## UQCR10 Rabbit Polyclonal Antibody

Catalog No: #54877

Package Size: #54877-1 50ul #54877-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	UQCR10 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UQCR10 (NP_037519.2).
Other Names	UQCR10;HSPC051;HSPC119;HSPC151;QCR9;UCCR7.2;UCRC
Accession No.	Uniprot:Q9UDW1GeneID:29796
Calculated MW	6kDa/7kDa
SDS-PAGE MW	7kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using UQCR10 antibody.

## Background

UCRC is a subunit of mitochondrial complex III (ubiquinol-cytochrome c reductase; EC 1.10.2.2), which forms the middle segment of the respiratory chain of the inner mitochondrial membrane (Schagger et al., 1995 [PubMed 8592474]).[supplied by OMIM, Mar 2008]

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.