

## OTUD7B Rabbit Polyclonal Antibody

Catalog No: #54883

Package Size: #54883-1 50ul #54883-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |   |
|-----------------------|---|
| Product Name          | OTUD7B Rabbit Polyclonal Antibody           |
| Host Species          | Rabbit                                      |
| Clonality             | Polyclonal                                  |
| Isotype               | IgG   |
| Purification          | Affinity purification                       |
| Applications          | WB,IHC                                      |
| Species Reactivity    | Human                                       |
| Immunogen Description | A synthetic Peptide of human OTUD7B         |
| Other Names           | OTUD7B;CEZANNE;ZA20D1                       |
| Accession No.         | Swiss Prot:Q6GQQ9GeneID:56957               |
| Calculated MW         | 88kDa/92kDa                                 |
| SDS-PAGE MW           | 93kDa                                       |
| Formulation           | Buffer: PBS with 0.02% sodium azide, pH7.3. |
| Storage               | Store at 4°C. Avoid freeze / thaw cycles.   |

## Application Details

WB □ 1:200 - 1:1000 IHC □ 1:20 - 1:50

## Images



Western blot analysis of extracts of MDA-MB-453 cells, using OTUD7B .

## Background

-

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.