

## TRIB3 Rabbit Polyclonal Antibody

Catalog No: #29679

Package Size: #29679-1 50ul #29679-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TRIB3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TRIB3 (NP_066981.2).
Other Names	TRIB3;C20orf97;NIPK;SINK;SKIP3;TRB3
Accession No.	Swiss Prot:Q96RU7GeneID:57761
Calculated MW	39kDa
SDS-PAGE MW	35kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using TRIB3 at 1:1000 dilution.

## Background

The protein encoded by this gene is a putative protein kinase that is induced by the transcription factor NF-kappaB. The encoded protein is a negative regulator of NF-kappaB and can also sensitize cells to TNF- and TRAIL-induced apoptosis. In addition, this protein can negatively regulate the cell survival serine-threonine kinase AKT1. Differential promoter usage and alternate splicing result in multiple transcript variants.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.