

NOLC1 Rabbit Polyclonal Antibody

Catalog No: #54936



Package Size: #54936-1 50ul #54936-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NOLC1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NOLC1 (NP_004732.2).
Other Names	NOLC1;NOPP130;NOPP140;NS5ATP13;P130
Accession No.	Swiss Prot:Q14978GeneID:9221
Calculated MW	73kDa/74kDa
SDS-PAGE MW	120kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:1000 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of U-2 OS cells using NOLC1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using NOLC1 at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.