

ENPEP Rabbit Polyclonal Antibody

Catalog No: #55034



Package Size: #55034-1 50ul #55034-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ENPEP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ENPEP (NP_001968.3).
Other Names	ENPEP;APA;CD249;gp160
Accession No.	Swiss Prot:Q07075GeneID:2028
Calculated MW	109kDa
SDS-PAGE MW	109-160kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200


Images




Western blot analysis of extracts of mouse lung, using ENPEP at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded rat kidney using ENPEP at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human prostate using ENPEP at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using ENPEP at dilution of 1:200 (40x lens).

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.