

TRIP10 Rabbit Polyclonal Antibody

Catalog No: #55092



Package Size: #55092-1 50ul #55092-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRIP10 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TRIP10 (NP_004231.1).
Other Names	TRIP10;CIP4;HSTP;STOT;STP;TRIP-10
Accession No.	Swiss Prot:Q15642GeneID:9322
Calculated MW	39kDa/52kDa/62kDa/67kDa/68kDa
SDS-PAGE MW	90kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details


WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

Images



Immunofluorescence analysis of U2OS cells using TRIP10 .
Blue: DAPI for nuclear staining.

Western blot analysis of extracts of various cell lines, using TRIP10 at 1:500 dilution.



Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.