

TACR3 Rabbit Polyclonal Antibody

Catalog No: #55148

Package Size: #55148-1 50ul #55148-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TACR3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TACR3 (NP_001050.1).
Other Names	TACR3;HH11;NK-3R;NK3R;NKR;TAC3RL
Accession No.	Uniprot:P29371GenelD:6870
Calculated MW	52kDa
SDS-PAGE MW	52kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TACR3 antibody.

Background

This gene belongs to a family of genes that function as receptors for tachykinins. Receptor affinities are specified by variations in the 5'-end of the sequence. The receptors belonging to this family are characterized by interactions with G proteins and 7 hydrophobic transmembrane regions. This gene encodes the receptor for the tachykinin neurokinin 3, also referred to as neurokinin B.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.