

# UBE2Z Rabbit Polyclonal Antibody

Catalog No: #55151



Package Size: #55151-1 50ul #55151-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	UBE2Z Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human UBE2Z (NP_075567.2).
Other Names	UBE2Z;HOYS7;USE1
Accession No.	Uniprot:Q9H832GeneID:65264
Calculated MW	28kDa/38kDa
SDS-PAGE MW	25kDa/38kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using UBE2Z antibody.

## Background

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.