

NGFRAP1 Rabbit Polyclonal Antibody

Catalog No: #55169



Package Size: #55169-1 50ul #55169-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NGFRAP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human NGFRAP1 (NP_996798.1).
Other Names	BEX3;Bex;DXS6984E;HGR74;NADE;NGFRAP1
Accession No.	Swiss Prot:Q00994GeneID:27018
Calculated MW	11kDa/12kDa
SDS-PAGE MW	24kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using NGFRAP1 at 1:1000 dilution.



Immunofluorescence analysis of MCF-7 cells using NGFRAP1

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.