

## TFEB Rabbit Polyclonal Antibody

Catalog No: #55175



Package Size: #55175-1 50ul #55175-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TFEB Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human TFEB.
Other Names	TFEB;ALPHATFEB;BHLHE35;TCFEB
Accession No.	Swiss Prot:P19484GeneID:7942
Calculated MW	43kDa/52kDa
SDS-PAGE MW	70kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:200


## Images



Immunofluorescence analysis of HeLa cells using TFEB at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using TFEB at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using TFEB at 1:1000 dilution.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.