

TCEAL8 Rabbit Polyclonal Antibody

Catalog No: #55345



Package Size: #55345-1 50ul #55345-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCEAL8 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TCEAL8 (NP_699164.1).
Other Names	TCEAL8;WEX3
Accession No.	Swiss Prot:Q8IYN2GeneID:90843
Calculated MW	13kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Images



Immunofluorescence analysis of A549 cells using TCEAL8 .
Blue: DAPI for nuclear staining.

Background

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternative splicing results in multiple transcript variants encoding a single isoform.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.