

USP26 Rabbit Polyclonal Antibody

Catalog No: #55396

Package Size: #55396-1 50ul #55396-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP26 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human USP26 (NP_114113.1).
Other Names	USP26
Accession No.	Uniprot:Q9BXU7GeneID:83844
Calculated MW	104kDa
SDS-PAGE MW	104kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using USP26 antibody.

Background

This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.