

DHX58 Rabbit Polyclonal Antibody

Catalog No: #55463



Package Size: #55463-1 50ul #55463-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DHX58 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DHX58 (NP_077024.2).
Other Names	DHX58;D11LGP2;D11lgp2e;LGP2;RLR-3
Accession No.	Swiss Prot:Q96C10GeneID:79132
Calculated MW	76kDa
SDS-PAGE MW	77kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using DHX58 at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using DHX58 at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.