

EIF3K Rabbit Polyclonal Antibody

Catalog No: #55642

Package Size: #55642-1 50ul #55642-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EIF3K Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human EIF3K (NP_037366.1).
Other Names	EIF3K;ARG134;EIF3-p28;EIF3S12;HSPC029;M9;MSTP001;PLAC-24;PLAC24;PRO1474;PTD001
Accession No.	Swiss Prot:Q9UBQ5GeneID:27335
Calculated MW	24kDa/25kDa
SDS-PAGE MW	25kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:100 - 1:200

Images



Western blot analysis of extracts of various cell lines, using EIF3K at 1:1000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.