

## CDC20 (phospho-Ser51) rabbit pAb

Catalog No: #13985



Package Size: #13985-1 50ul #13985-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CDC20 (phospho-Ser51) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	This antibody detects endogenous levels of Human Mouse CDC20 (phospho-Ser51)
Immunogen Description	Synthesized phosho peptide around human CDC20 (Ser51)
Other Names	Cell division cycle protein 20 homolog (p55CDC)
Accession No.	Swiss Prot:Q12834GeneID:991
SDS-PAGE MW	50
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

## Background

cell division cycle 20(CDC20) Homo sapiens CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation. [provided by RefSeq, Jul 2008],

Note: This product is for in vitro research use only and is not intended for use in humans or animals.