

# Recombinant Human Fibronectin type III domain-containing protein 3B(FNDC3B) ,partial

Catalog No: #AP0103

Package Size: #AP0103-1 20ug #AP0103-2 100ug #AP0103-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Fibronectin type III domain-containing protein 3B(FNDC3B) ,partial
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Name	His-tagged
Other Names	FNDC3B; FAD104; NS5ABP37; UNQ2421/PRO4979/PRO34274Fibronectin type III domain-containing protein 3B; Factor for adipocyte differentiation 104; HCV NS5A-binding protein 37
Tag Info	N-terminal His-tagged
Target Sequence	278-570aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

Greater than 85% as determined by SDS-PAGE.

## Background

May be a positive regulator of adipogenesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.