

# Recombinant Human Ras-related protein Rab-33B(RAB33B)

Catalog No: #AP0105

Package Size: #AP0105-1 20ug #AP0105-2 100ug #AP0105-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Ras-related protein Rab-33B(RAB33B)
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Name	His-tagged
Other Names	DKFZP434G099; MGC138182; RAB 33B; RAB33B; RAB33B member RAS oncogene family; Ras related protein Rab 33B; Ras related protein Rab33B; Ras-associated protein RAB33B; Ras-related protein Rab-33B; RB33B_HUMAN; SMC2
Tag Info	N-terminal His-tagged
Target Sequence	1-229aa
Concentration	0.1-1.0 mg/mL, by the Bradford Method.
Formulation	Tris-based buffer, 50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

Greater than 85% as determined by SDS-PAGE.

## Background

Protein transport. Acts, in coordination with RAB6A, to regulate intra-Golgi retrograde trafficking. It is involved in autophagy, acting as a modulator of autophagosome formation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.