

Recombinant Human Protein FAM190B(FAM190B) ,partial

Catalog No: #AP0110



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP0110-1 20ug #AP0110-2 100ug #AP0110-3 1mg

Description

Product Name	Recombinant Human Protein FAM190B(FAM190B) ,partial
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Sequence	231-328aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.