

Human PSEN1 / Presenilin-1 ELISA Kit

Catalog No: #EK21005



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human PSEN1 / Presenilin-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	PS-1,Protein S182,AD3,PS1,PSNL1
Accession No.	Uniprot:P49768
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:21pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.