

# Human VASH2 / Tubuliny-Tyr carboxypeptidase 2 ELISA Kit

Catalog No: #EK21006



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human VASH2 / Tubuliny-Tyr carboxypeptidase 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	VASHL, Vasohibin-2, Vasohibin-like protein
Accession No.	Uniprot:Q86V25
Storage	6 months

## Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.24ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.