

Human FDXR / NADPH:adrenodoxin oxidoreductase, mitochondrial ELISA Kit

Catalog No: #EK21011



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human FDXR / NADPH:adrenodoxin oxidoreductase, mitochondrial ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	AR, Ferredoxin--NADP(+) reductase, Ferredoxin reductase, ADXR
Accession No.	Uniprot:P22570
Storage	6 months

Application Details

Detection Range:0.625-40ng/mL Assay Type:Sandwich Sensitivity:0.284ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.