

Rat Hfe / Hereditary hemochromatosis protein homolog ELISA Kit

Catalog No: #EK21012



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Hfe / Hereditary hemochromatosis protein homolog ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	RT
Other Names	RT1-CAFE
Accession No.	Uniprot:O35799
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.27ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.