

Mouse Spam1 / Hyaluronidase PH-20 ELISA Kit

Catalog No: #EK21015



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Spam1 / Hyaluronidase PH-20 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	Hyal-PH20, Hyaluronoglucosaminidase PH-20, Sperm adhesion molecule 1, Sperm surface protein PH-20, Ph20, Spam
Accession No.	Uniprot:P48794
Storage	6 months

Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.064ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.