

Human CNTNAP4 / Contactin-associated protein-like 4 ELISA Kit

Catalog No: #EK21018



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CNTNAP4 / Contactin-associated protein-like 4 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cell recognition molecule Caspr4,CASPR4,KIAA1763
Accession No.	Uniprot:Q9C0A0
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.056ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.