

# Mouse Mtap / S-methyl-5'-thioadenosine phosphorylase ELISA Kit

Catalog No: #EK21019



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Mtap / S-methyl-5'-thioadenosine phosphorylase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	5'-methylthioadenosine phosphorylase,MTA phosphorylase
Accession No.	Uniprot:Q9CQ65
Storage	6 months

## Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.