

Human CSRP3 / Cysteine and glycine-rich protein 3 ELISA Kit

Catalog No: #EK21020



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human CSRP3 / Cysteine and glycine-rich protein 3 ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Human |
| Other Names | Cardiac LIM protein, Cysteine-rich protein 3, CRP3, LIM domain protein, cardiac, Muscle LIM protein, CLP, MLP |
| Accession No. | Uniprot:P50461 |
| Storage | 6 months |

Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.078ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.