

# Human METTL18 / Histidine protein methyltransferase 1 homolog ELISA Kit

Catalog No: #EK21021



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human METTL18 / Histidine protein methyltransferase 1 homolog ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Arsenic-transactivated protein 2,AsTP2,Methyltransferase-like protein 18,ASTP2,C1orf156
Accession No.	Uniprot:O95568
Storage	6 months

## Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:36pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.