

Human WNT3 / Proto-oncogene Wnt-3 ELISA Kit

Catalog No: #EK21022



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human WNT3 / Proto-oncogene Wnt-3 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Proto-oncogene Int-4 homolog,INT4
Accession No.	Uniprot:P56703
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.14ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.