

# Human CFHR5 / Complement factor H-related protein 5 ELISA Kit

Catalog No: #EK21026



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CFHR5 / Complement factor H-related protein 5 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	FHR-5,CFHL5,FHR5
Accession No.	Uniprot:Q9BXR6
Storage	6 months

## Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:30pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.