

Human MSMP / Prostate-associated microseminoprotein ELISA Kit

Catalog No: #EK21027



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human MSMP / Prostate-associated microseminoprotein ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	PC3-secreted microprotein,PSMP
Accession No.	Uniprot:Q1L6U9
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.085ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.