

Human HNRNPH2 / Heterogeneous nuclear ribonucleoprotein H2 ELISA Kit

Catalog No: #EK21028



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human HNRNPH2 / Heterogeneous nuclear ribonucleoprotein H2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	hnRNP H2,FTP-3,Heterogeneous nuclear ribonucleoprotein H',hnRNP H',FTP3,HNRPH2
Accession No.	Uniprot:P55795
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.084ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.