

Human LRRC10B / Leucine-rich repeat-containing protein 10B ELISA Kit

Catalog No: #EK21058



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human LRRC10B / Leucine-rich repeat-containing protein 10B ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:A6NIK2
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.158ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.