

Human NEU1 / Sialidase-1 ELISA Kit

Catalog No: #EK21059



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human NEU1 / Sialidase-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Acetylneuraminyl hydrolase,G9 sialidase,Lysosomal sialidase,N-acetyl-alpha-neuraminidase 1,NANH
Accession No.	Uniprot:Q99519
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.22ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.