

Human BNIP3L / BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like ELISA Kit

Catalog No: #EK21060



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human BNIP3L / BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Adenovirus E1B19K-binding protein B5,BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A,NIP3-like protein X,NIP3L,BNIP3A,BNIP3H,NIX
Accession No.	Uniprot:O60238
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.14ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.