

Human CBL / E3 ubiquitin-protein ligase CBL ELISA Kit

Catalog No: #EK21063



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CBL / E3 ubiquitin-protein ligase CBL ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Casitas B-lineage lymphoma proto-oncogene,Proto-oncogene c-Cbl,RING finger protein 55,RING-type E3 ubiquitin transferase CBL,Signal transduction protein CBL,CBL2,RNF55
Accession No.	Uniprot:P22681
Storage	6 months

Application Details

Detection Range:1.56-100ng/mL Assay Type:Sandwich Sensitivity:0.81ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.