

Human LYN / Tyrosine-protein kinase Lyn ELISA Kit

Catalog No: #EK21064



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human LYN / Tyrosine-protein kinase Lyn ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Lck/Yes-related novel protein tyrosine kinase,V-yes-1 Yamaguchi sarcoma viral related oncogene homolog,p53Lyn,p56Lyn,JTK8
Accession No.	Uniprot:P07948
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.057ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.