

Human ASPG / 60 kDa lysophospholipase ELISA Kit

Catalog No: #EK21066



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ASPG / 60 kDa lysophospholipase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	C14orf76,Lysophospholipase-transacylase
Accession No.	Uniprot:Q86U10
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.039ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.