

Human FOLH1 / Glutamate carboxypeptidase 2 ELISA Kit

Catalog No: #EK21067



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human FOLH1 / Glutamate carboxypeptidase 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cell growth-inhibiting gene 27 protein, Folate hydrolase 1, Folylpoly-gamma-glutamate carboxypeptidase, FGCP, Glutamate carboxypeptidase II, GCPII, Membrane glutamate carboxypeptidase, mGCP, N-acetylated-alpha-linked acidic dipeptidase I, NAALADase I, Prostate-specific membrane antigen, PSM, Pteroylpoly-gamma-glutamate carboxypeptidase, FOLH, NAALAD1, PSM, PSMA, GIG27
Accession No.	Uniprot:Q04609
Storage	6 months

Application Details

Detection Range: 0.625-40ng/mL Assay Type: Sandwich Sensitivity: 0.295ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.