

Human ZC3H8 / Zinc finger CCCH domain-containing protein 8 ELISA Kit

Catalog No: #EK21068



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ZC3H8 / Zinc finger CCCH domain-containing protein 8 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	ZC3HDC8
Accession No.	Uniprot:Q8N5P1
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.057ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.