

Human NRIP1 / Nuclear receptor-interacting protein 1 ELISA Kit

Catalog No: #EK21069



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human NRIP1 / Nuclear receptor-interacting protein 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Nuclear factor RIP140, Receptor-interacting protein 140
Accession No.	Uniprot:P48552
Storage	6 months

Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.13ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.