

Human MFAP3 / Microfibril-associated glycoprotein 3 ELISA Kit

Catalog No: #EK21071



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human MFAP3 / Microfibril-associated glycoprotein 3 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:P55082
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:36pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.