

Human CCAR2 / Cell cycle and apoptosis regulator protein 2 ELISA Kit

Catalog No: #EK21072



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CCAR2 / Cell cycle and apoptosis regulator protein 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cell division cycle and apoptosis regulator protein 2, DBIRD complex subunit KIAA1967, Deleted in breast cancer gene 1 protein, DBC-1, NET35, p30 DBC, DBC1, KIAA1967
Accession No.	Uniprot:Q8N163
Storage	6 months

Application Details

Detection Range: 0.78-50ng/mL Assay Type: Sandwich Sensitivity: 0.33ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.