

Human TBXA2R / Thromboxane A2 receptor ELISA Kit

Catalog No: #EK21073



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human TBXA2R / Thromboxane A2 receptor ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	TXA2-R, Prostanoid TP receptor
Accession No.	Uniprot:P21731
Storage	6 months

Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.039ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.