

# Human MT-ND1 / NADH-ubiquinone oxidoreductase chain 1 ELISA Kit

Catalog No: #EK21074



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human MT-ND1 / NADH-ubiquinone oxidoreductase chain 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	NADH dehydrogenase subunit 1,MTND1,NADH1,ND1
Accession No.	Uniprot:P03886
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.144ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.