

Mouse Lingo1 / Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1 ELISA Kit

Catalog No: #EK21076

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Lingo1 / Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	Leucine-rich repeat neuronal protein 1, Leucine-rich repeat neuronal protein 6A, Lern1, Lrrn6a
Accession No.	Uniprot:Q9D1T0
Storage	6 months

Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.11ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.