

# Human ECD / Protein ecdysoneless homolog ELISA Kit

Catalog No: #EK21080



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ECD / Protein ecdysoneless homolog ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Human suppressor of GCR two,hSGT1
Accession No.	Uniprot:O95905
Storage	6 months

## Application Details

Detection Range:31.2-2000pg/mL Assay Type:Sandwich Sensitivity:12.5 pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.