

Human MPP1 / 55 kDa erythrocyte membrane protein ELISA Kit

Catalog No: #EK21081



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human MPP1 / 55 kDa erythrocyte membrane protein ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	p55,Membrane protein, palmitoylated 1,DXS552E,EMP55
Accession No.	Uniprot:Q00013
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.36ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.