

# Human SNRPB / Small nuclear ribonucleoprotein-associated proteins B and B' ELISA Kit



Catalog No: #EK21082

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Human SNRPB / Small nuclear ribonucleoprotein-associated proteins B and B' ELISA Kit |
| Brief Description  | ELISA Kit  |
| Applications       | ELISA  |
| Species Reactivity | Human  |
| Other Names        | snRNP-B, Sm protein B/B', Sm-B/B', COD, SNRPB1                                       |
| Accession No.      | Uniprot:P14678   |
| Storage            | 6 months   |

## Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.051ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.