

# Human BRD2 / Bromodomain-containing protein 2 ELISA Kit

Catalog No: #EK21083



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human BRD2 / Bromodomain-containing protein 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Really interesting new gene 3 protein, KIAA9001, RING3, O27.1.1
Accession No.	Uniprot:P25440
Storage	6 months

## Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.11ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.